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EN

## **Warnings and safety information**

All information regarding the use of the product can be found in the instruction manual. Before using it, read its contents and follow the instructions contained therein.

**Also read the following information before use:**

### **Warnings for use**

#### **Risk of product damage:**

- Use the devices only as intended, avoiding use on the body in a manner inconsistent with the instructions.
- Non-waterproof devices should not be immersed in water or used in high humidity conditions.
- Do not expose products to extreme temperatures or leave them in the sun for long periods of time.
- For battery-powered devices, make sure to charge only with the supplied cable or compatible adapter.

#### **User Hazard:**

- Avoid using devices on irritated, damaged skin, open wounds, active inflammation or allergic reactions.
- Do not use devices that emit microcurrents, ultrasound or heat if you have a pacemaker, metal implants or are pregnant.
- If you experience sensitization, discomfort or other worrisome symptoms, stop using the device immediately and consult a dermatologist.

## Information on proper use

### Recommendations for use:

- Before using devices such as massagers or mustache vacuums, cleanse your face thoroughly and apply appropriate cosmetics (such as gel or water) to support the treatment.
- Deep-cleansing devices, such as cavitation peels, should be used a maximum of 2-3 times a week to avoid excessive skin irritation.
- Use the devices according to the manufacturer's recommendations, not exceeding the recommended duration of action on one area of the face (usually 1-2 minutes).
- Choose the appropriate intensity level, especially if you have sensitive skin.

### Maintenance and storage:

- After each use, clean the device with a damp cloth or other tool according to the manufacturer's recommendations. For waterproof devices, you can rinse them under running water.
- Regularly check the condition of accessories, such as scrub or brush tips, and replace them according to the manufacturer's recommendations.
- Store devices in a dry, clean place, avoiding contact with moisture and dust.

## Additional precautions

### Customization to the user:

- Choose the right device for your skin type and individual needs. Functions such as heating or micro-current massage may not be advisable for delicate or vascular skin.
- Use of devices by people with dermatological problems, such as rosacea, should be done in consultation with a doctor.

### Damage and service:

- If you notice any damage to the device (e.g., cracks in the case, problems with switching on, unstable operation), immediately stop using and contact the manufacturer's authorized service center.
- Use only original components supplied by the manufacturer to replace parts and accessories.
- Regularly check the performance of rechargeable devices and charge them as recommended to avoid performance degradation or risk of damage.

## Statement of compliance

The products have been designed in accordance with the safety requirements of Regulation (EU) 2023/988. They meet the relevant technical standards, including those for safety of use, electromagnetic compatibility and durability of materials.

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