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EN

## **Warnings and safety information**

All information regarding the use of the product can be found in the instruction manual. Before using it, read its contents and follow the instructions contained therein.

**Also read the following information before use:**

### **Warnings for use**

#### **Risk of product damage:**

- Use the irrigator as intended and as recommended by the manufacturer.
- Avoid using water at a higher temperature than recommended in the instructions to prevent damage to internal components.
- Do not immerse the device in water unless it is marked as fully waterproof.
- Regularly check the condition of terminals and power cords - worn or damaged components should be replaced immediately.

#### **Risks to the user:**

- Irrigators are not for people with open mouth wounds, gum infections or sensitive teeth, unless recommended by a dentist.
- The devices should not be used by children without adult supervision.
- Improper use, such as directing the water jet directly at sensitive tissues, can lead to irritation.

## **Information on proper use**

#### **Recommendations for use:**

- Choose an irrigator tip appropriate for your needs - most devices come with a variety of tips that can be used for interdental cleaning, gum care or other specific applications, according to the manufacturer's recommendations.
- Set the appropriate level of water pressure, starting at the lowest level to avoid discomfort.
- Wash the irrigator tips regularly and replace them according to the manufacturer's recommendations (usually every 3-6 months).
- Remove water from the device after each use to prevent bacteria buildup.

#### **Maintenance and storage:**

- Clean the irrigator as recommended in the instruction manual. Avoid using aggressive detergents.
- Store the device in a dry place, away from sources of moisture and extreme temperatures.

- If the irrigator has a charging function, use only the original cable provided by the manufacturer.

## **Additional precautions**

### **Adaptation to the user:**

- People with braces should consult their dentist before using the irrigator to determine appropriate use.
- Children and the elderly should only use the device under supervision.

### **Damage and service:**

- If you notice any damage to the device (e.g., leaks, water jet problems), discontinue use and contact the manufacturer's service department.
- Use only original components approved by the manufacturer for repairs and replacement of parts.

## **Statement of conformity**

The products have been designed in accordance with the safety requirements of Regulation (EU) 2023/988. They meet the relevant technical standards, including safety of use, durability of materials and protection of consumer health.

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