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EN

Warnings and safety information

All information regarding the use of the product can be found in the instruction manual. Before using it, read its contents and follow the instructions contained therein.

Also read the following information before use:

Warnings for use

Risk of product damage:

- Use sonic toothbrushes only as intended and as recommended by the manufacturer.
- Regularly check the condition of the device, including the tips and charger. In case of visible damage, such as cracks or abrasions, stop using immediately.
- Do not completely immerse the devices in water unless they are marked waterproof.
- Store the device in a dry place to avoid damage from excessive moisture.

User Hazard:

- Avoid using sonic toothbrushes on irritated or damaged oral skin.
- Do not direct the irrigator's stream directly onto injured or sensitive areas of the mouth.
- People with braces or implants should consult a dentist before using a sonic toothbrush.
- Children should only use sonic toothbrushes under adult supervision.

Risks associated with electronic components:

- Use only chargers and cables supplied by the manufacturer to avoid the risk of overheating batteries.
- Do not use the device if it shows signs of electrical damage (e.g. erratic operation, heating up).
- Connect USB-charged products only to compatible power sources.

Information on proper use**Recommendations for use:**

- Use tips appropriate to your needs.
- Always use the right toothpaste, preferably recommended by the toothbrush manufacturer.
- Set the toothbrush mode to gentle if you are just starting to use sonic devices.
- Adjust the strength of the irrigator's water jet to suit your needs. Choose the gentle mode to start with, especially if you have sensitive gums.
- Regularly refill your irrigator's reservoir with water or a specialized oral hygiene liquid.
- For sonic body brushes, use the appropriate pads for your chosen area of care.

Maintenance and storage:

- Rinse your toothbrush thoroughly after each use and dry it before storing.
- Replace your toothbrush tips regularly - this is recommended every 3 months or more often for heavy use.
- For devices with built-in batteries, charge them according to the manufacturer's recommendations to avoid damage.

Additional precautions**Adaptation to the user:**

- People with dental problems, such as tooth sensitivity or gum disease, should consult a dentist before using sonic toothbrushes or irrigators.
- Regularly observe the reactions of your skin and gums. If irritation occurs, discontinue use of the device and consult your doctor.
- Use devices according to their intended use - sonic body brushes should not be used on facial skin, and irrigators should only be used in the mouth.

Damage and service:

- If you notice any problems with the operation of the device (e.g., loss of power, leaking irrigator, damaged tips), stop using it immediately and contact an authorized service center.
- Use only original products recommended by the manufacturer to replace tips, accessories or spare parts.
- Do not attempt to repair the devices yourself - repairs should be performed only by specialists at authorized service centers.

Statement of conformity

The products have been designed in accordance with the safety requirements of Regulation (EU) 2023/988. They meet the relevant technical standards, including those for safety of use, electromagnetic compatibility and durability of materials.

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